

SACROSANCT International

BONE FORCEPS



www.sacrointl.com



PROFILE

SACROSANCT INTERNATIONAL is a leading, family owned company, dedicated to gives you the personalized service you deserve. We have been producing & supplying the finest Health Care Instruments all over the world. Our Production Houses are located in Sialkot, The City of Health Care Instrument's Manufacturers.

Sacrosanct International offers a broad line of General Surgical Instruments, T.C Instruments, Dental Instruments, Eye Instruments & Student Dissecting Sets.

This Catalogue is just for introduction , if you not find the items of your interest , send us the images & description of your interested items, We would be Happy to produce for you the items of your interest. Please feel free if you have any question.

Thanks & Best Regards
Slah-ud-Din
Proprietor



BONE FORCEPS



SI-BF-001

Vogen Sesamoid Clamp

Designed to hold a sesamoid during dissection. 4"



SI-BF-002

Stagbeetle Forceps 4 3/4"

Maximum opening with ratchet engaged: 14 mm



SI-BF-003

Plate and Bone Holding Forceps 5"

One pointed tip, one footplate



SI-BF-004

Bone Reduction Forceps 5"

Curved pointed tips



SI-BF-005

Bone Reduction Forceps 5 1/4"

Curved pointed tips



SI-BF-006

Curved Serrated jaws 6"

Curved pointed tips



SI-BF-007

Forceps Curved with 5 1/4"

Curved one pointed tip
one step-pointed tip



SI-BF-008

Bone Reduction Forceps 6 3/4"

Curved Serrated jaws



SI-BF-009

Locke Phalangeal 6"
Serrated jaws with 1x2 teeth





BONE FORCEPS



SI-BF-0010
Reill Reduction
Forceps 7" Curved jaw



SI-BF-0011
Bone Reduction
Forceps 7" One step-pointed
tip 2" double ratchet



SI-BF-0012
Bone Reduction
Forceps 8"
2" double ratchet



SI-BF-0013
Bone Reduction
Forceps 8"
With speed lock Max



SI-BF-0014
Lewin Bone Forceps 7"
Slightly curved overlapping
serrated jaws



SI-BF-0015
Malleolar Forceps 8"
With speed lock Max



SI-BF-0016
Dingman Forceps 7 1/2"
Delicate serrated jaws for
small bone fragments



SI-BF-0017
Sesamoidectomy
Clamp 6 1/2"



SI-BF-0018
Fibula Forceps 8 1/2"
With speedlock angled tips





BONE FORCEPS



SI-BF-0019
Patella Forceps 7"
Speedlock 2x2 sharp teeth



SI-BF-0020
Bone Holding Forceps 7"
With speedlock. Curved serrated
jaws



SI-BF-0021
Plate Holding Forceps 5 1/2"



SI-BF-0022
Plate Holding Forceps 7 1/2"
For positioning 2.7mm and 3.5mm
plates



SI-BF-0023
Verbrugge Forceps 6 1/2"
Reversed jaw with long ratchet



SI-BF-0024
Verbrugge Forceps 7"
With short ratchet



SI-BF-0025
Semb Bone Forceps 7 1/2"
Angled on side Deep teeth,
without ratchet



SI-BF-0026
Semb Bone Forceps 7 1/2"
Angled on side Deep teeth,
with ratchet



SI-BF-0027
Van Buren Bone Forceps 8 1/2"
Angled on side Serrated jaws





BONE FORCEPS



SI-BF-0028
Fergusson Bone Forceps 8"
Deep 2x2 teeth



SI-BF-0029
Langenbeck Bone Forceps 8"
Partially serrated. 2x2 teeth



SI-BF-0030
Kern Bone Forceps
5, 1/2", 6, 1/2", 9, 1/2"
Serrated jaws. 2x2 teeth



SI-BF-0031
Lane Bone Forceps
13" light, 13" standard, 18" heavy



SI-BF-0032
Lane Bone Forceps
13" light, 13" standard, 18" heavy



SI-BF-0033
Lowman Bone
Clamp 1x1 jaws.
4, 3/4", 7, 1/4" 8"



SI-BF-0034
Lowman Bone
Clamp 1x2 jaws.
5", 7, 1/4" 8"



SI-BF-0035
Lambert-Lowman Bone
Clamp 2x2 jaws.
7, 1/4" 8"



SI-BF-0036
Ulrich Bone Forceps
Straight. with speedlock
7, 1/4", 9", 9 1/2" 11"





BONE FORCEPS



SI-BF-0037
Farabeuf Forceps
Adjustable jaw 9", 10"



SI-BF-0038
Farabeuf Lambotte Forceps
Adjustable jaw with ratchet 7 1/2"



SI-BF-0039
Lambotte Bone Forceps
Adjustable jaw with swivel
head, with ratchet
8, 1/2, 10, 1/2" 11, 1/2"



SACROSANCT International

📞 +92 333 8682068

✉ info@sacrointl.com

🌐 www.sacrointl.com

